



LOWER EXTREMITY TRAUMA PLATES



For our complete range of systems, please tab or scan the QR code:





LOWER EXTREMITY
TRAUMA PLATES

working for patient comfort



Stabilize where force meets anatomy.



ortho

LOWER EXTREMITY
TRAUMA PLATES

For foot plate systems refer to our Mini Fragment Plating System Catalog and also for our complete range of systems, please [visit the link](#) - tab or scan the QR code:



Ankle



Tibia



Femur & Knee



Hip & Pelvic



ANKLE



Locking Calcaneal Plate




Plate Thickness	2.0mm
Plate Shaft Width	
Plate Head Width	
Shaft Hole Spacing	
Material	Titanium

Compatible with:
3.5 mm Locking Screw
3.5 mm Cortical Screw
4.0mm Cancellous Screw

Ref. No	Hole	Direction	Plate Length
P00482	Small	Left	57.90 mm
P00483	Medium	Left	72.80 mm
P00115	Large	Left	77 mm
P00484	Small	Right	57.90 mm
P00485	Medium	Right	72.80 mm
P00116	Large	Right	77 mm

Distal Anterior-Lateral Tibia Locking Plate



Plate Thickness	3.50mm
Plate Shaft Width	13mm
Plate Head Width	32mm
Shaft Hole Spacing	16mm
Material	Titanium

Compatible with:
3.5mm Locking Screw
3.5mm Cortical Screw
4.0mm Cancellous Screw
3.5mm Spacer

Ref. No	Hole	Direction	Plate Length
P00664	4	Left	84.3 mm
P00665	4	Right	84.3 mm
P00666	6	Left	116.3 mm
P00667	6	Right	116.3 mm
P00668	8	Left	148.3 mm
P00669	8	Right	148.3 mm
P00670	10	Left	180.3 mm
P00671	10	Right	180.3 mm
P00672	12	Left	212.3 mm
P00673	12	Right	212.3 mm
P00674	14	Left	228.3 mm
P00675	14	Right	228.3 mm
P00676	16	Right	260.3 mm
P00677	16	Left	260.3 mm
P00678	18	Right	292.3 mm
P00679	18	Left	292.3 mm

Locking LCP Distal Lateral Fibula Plate



Plate Thickness	2.50mm
Plate Shaft Width	9mm
Plate Head Width	13.50mm
Shaft Hole Spacing	13mm
Material	Titanium

Compatible with:
2.4mm/2.7mm/3.5mm Locking Screw
3.5mm Cortical Screw
4.0mm Cancellous Screw
3.5mm Spacer

Ref. No	Hole	Direction	Plate Length
P00229	3	Left	73.0 mm
P00230	4	Left	86.0 mm
P00231	5	Left	99.0 mm
P00232	6	Left	112.0 mm
P00233	7	Left	125.0 mm
P00571	8	Left	138.0 mm
P00572	9	Left	151.0 mm
P00234	3	Right	73.0 mm
P00235	4	Right	86.0 mm
P00236	5	Right	99.0 mm
P00237	6	Right	112.0 mm
P00238	7	Right	125.0 mm
P00573	8	Right	138.0 mm
P00574	9	Right	151.0 mm

Locking Distal Fibula Plate 3.5mm




Plate Thickness	2.0mm
Plate Shaft Width	9.60mm
Plate Head Width	14.50mm
Shaft Hole Spacing	13mm
Material	Titanium

Compatible with:
3.5mm Locking Screw
3.5mm Cortical Screw
4.0mm Cancellous Screw
3.5mm Spacer

Ref. No	Hole	Direction	Plate Length
P00239	2	Left	60.4 mm
P00240	4	Left	86.4 mm
P00241	6	Left	112.4 mm
P00242	8	Left	138.4 mm
P00243	10	Left	164.4 mm
P00244	2	Right	60.4 mm
P00245	4	Right	86.4 mm
P00246	6	Right	112.4 mm
P00247	8	Right	138.4 mm
P00248	10	Right	164.4 mm

3.5 mm Locking LCP Distal Medial Tibia Plate



Compatible with:
3.5mm Locking Screw
3.5mm Cortical Screw
4.0mm Cancellous Screw
3.5 mm Spacer

Material Titanium

Locking Distal Posterolateral Fibula Plate



Compatible with:
2.4mm/2.7mm/3.5mm Locking Screw
3.5mm Cortical Screw
4.0mm Cancellous Screw

Material Titanium

Ref. No	Hole	Direction
P00508	4	Left
P00509	6	Left
P00510	8	Left
P00511	10	Left
P00512	12	Left
P00513	14	Left
P00514	16	Left
P00515	18	Left
P00516	20	Left
P00517	4	Right
P00518	6	Right
P00519	8	Right
P00520	10	Right
P00521	12	Right
P00522	14	Right
P00523	16	Right
P00524	18	Right
P00525	20	Right

Ref. No	Hole	Direction
P00575	3	Left
P00576	4	Left
P00577	5	Left
P00578	6	Left
P00579	7	Left
P00580	8	Left
P00581	9	Left
P00582	3	Right
P00583	4	Right
P00584	5	Right
P00585	6	Right
P00586	7	Right
P00587	8	Right
P00588	9	Right

Distal Posterior Tibia T-Plate



Compatible with:
3.5mm Locking Screw
3.5mm Cortical Screw
4.0mm Cancellous Screw

Material Titanium

Distal Posterior Tibia L-Plate



Compatible with:
3.5mm Locking Screw
3.5mm Cortical Screw
4.0mm Cancellous Screw

Material Titanium

Ref. No	Hole
P01200	3
P01201	4
P01202	5

Ref. No	Hole	Direction
P01194	3	Right
P01195	4	Right
P01196	5	Right
P01197	3	Left
P01198	4	Left
P01199	5	Left

Locking LCP Distal Lateral Fibula Plate



Compatible with:
3.5mm Locking Screw
3.5mm Cortical Screw
4.0mm Cancellous Screw
3.5 mm Spacer

Material Titanium

Ref. No	Hole	Direction
P00854	3	Universal
P00855	4	Universal
P00856	5	Universal
P00857	6	Universal
P00858	7	Universal
P00859	9	Universal

3.5 mm Locking Distal Fibula Plate



Compatible with:
3.5mm Locking Screw
3.5mm Cortical Screw
4.0mm Cancellous Screw
3.5 mm Spacer

Material Titanium

Ref. No	Hole	Direction
P00860	2	Universal
P00861	4	Universal
P00862	6	Universal
P00863	8	Universal
P00864	10	Universal

TIBIA



Locking Proximal Tibia Osteotomy Plate




Plate Thickness	3.50mm
Plate Shaft Width	15.90mm
Plate Head Width	33.80mm
Shaft Hole Spacing	
Material	Titanium

Compatible with:
5.0mm Locking Screw
4.5mm Cortical Screw
5.0mm Spacer

Ref. No	Hole	Plate Length
P00680	4	114.3 mm

4.5 mm Locking Proximal Medial Tibia T-Plate




Plate Thickness	2.50mm
Plate Shaft Width	15mm
Plate Head Width	34mm
Shaft Hole Spacing	16mm
Material	Titanium

Compatible with:
5.0mm Locking Screw
4.5mm Cortical Screw
5.0mm Spacer

Ref. No	Hole	Plate Length
P00698	4	105 mm
P00699	6	137 mm
P00700	8	169 mm

Locking LCP Proximal Medial Tibia Plate




Plate Thickness	2.50mm
Plate Shaft Width	15mm
Plate Head Width	33mm
Shaft Hole Spacing	16mm
Material	Titanium

Compatible with:
5.0mm Locking Screw
4.5mm Cortical Screw
5.0mm Spacer

Ref. No	Hole	Direction
P00546	3	Left
P00223	4	Left
P00547	5	Left
P00224	6	Left
P00548	7	Left
P00225	8	Left
P00549	3	Right
P00226	4	Right
P00550	5	Right
P00227	6	Right
P00551	7	Right
P00228	8	Right

Locking Proximal Lateral Tibia Plate



Material: Titanium

Compatible with:
3.5mm Locking Screw
3.5mm Cortical Screw
4.0mm Cancellous Screw
3.5mm Spacer

Ref. No	Hole	Direction
P00528	3	Left
P00529	5	Left
P00530	7	Left
P00531	9	Left
P00532	11	Left
P00533	13	Left
P00534	15	Left
P00535	17	Left
P00536	3	Right
P00537	5	Right
P00538	7	Right
P00539	9	Right
P00540	11	Right
P00541	13	Right
P00542	15	Right
P00543	17	Right

Tibia

3.5 mm Locking Proximal Medial Tibia Plate



Compatible with:
3.5mm Locking Screw
3.5mm Cortical Screw
4.0mm Cancellous Screw
3.5mm Spacer

Material Titanium

Ref. No	Hole	Direction
P00552	3	Left
P00553	5	Left
P00554	7	Left
P00555	9	Left
P00556	11	Left
P00557	13	Left
P00558	15	Left
P00559	3	Right
P00560	5	Right
P00561	7	Right
P00562	9	Right
P00563	11	Right
P00564	13	Right
P00565	15	Right

3.5 mm Locking Posteromedial Tibia Plate



Compatible with:
3.5mm Locking Screw
3.5mm Cortical Screw
4.0mm Cancellous Screw
3.5mm Spacer

Material Titanium

Ref. No	Hole	Direction
P00566	3	
P00567	5	
P00568	7	
P00569	9	
P00570	11	

Locking Proximal Tibia Plate




Plate Thickness	5.0mm
Plate Shaft Width	15.50mm
Plate Head Width	22mm
Shaft Hole Spacing	17mm
Material	Titanium

Compatible with:
5.0mm Locking Screw
4.5 mm Cortical Screw
5.0mm Cancellous Screw
5.0mm Locking Cancellous Screw
5.0 mm Spacer

Ref. No	Hole	Direction	Plate Length
P00199	3	Left	104.50 mm
P00200	5	Left	138.50 mm
P00201	7	Left	172.5 mm
P00202	9	Left	206.5 mm
P00203	11	Left	240.5 mm
P00204	13	Left	274.5 mm
P00526	15	Left	309.5 mm
P00205	3	Right	104.5 mm
P00206	5	Right	138.5 mm
P00207	7	Right	172.5 mm
P00208	9	Right	206.5 mm
P00209	11	Right	240.5 mm
P00210	13	Right	274.5 mm
P00527	15	Right	309.5 mm

4.5mm Locking Proximal Medial Tibia T-Plate



Plate Thickness	2.50mm
Plate Shaft Width	16mm
Plate Head Width	38mm
Shaft Hole Spacing	16mm
Material	Titanium

Compatible with:
5.0mm Locking Screw
4.5 mm Cortical Screw
5.0mm Cancellous Screw
5.0mm Locking Cancellous Screw
5.0 mm Spacer

Ref. No	Hole	Direction	Plate Length
P00179	4	Left	96.52 mm
P00180	6	Left	124.52 mm
P00181	8	Left	156.52 mm
P00182	10	Left	188.52 mm
P00183	4	Right	96.52 mm
P00184	6	Right	124.52 mm
P00185	8	Right	156.52 mm
P00186	10	Right	188.52 mm

Locking 4.5 Narrow Compression Plate



Plate Thickness	5.0mm
Plate Shaft Width	13.80mm
Plate Head Width	
Shaft Hole Spacing	18mm
Material	Titanium

Compatible with:
5.0mm Locking Screw
4.5 mm Cortical Screw
5.0mm Cancellous Screw

Ref. No	Hole	Plate Length
P00276	4	84 mm
P00277	5	102 mm
P00278	6	120 mm
P00279	7	138 mm
P00280	8	156 mm
P00281	9	174 mm
P00282	10	192 mm
P00283	11	210 mm
P00284	12	228 mm
P00285	13	246 mm
P00286	14	264 mm
P00287	15	282 mm
P00288	16	300 mm

FEMUR & KNEE



Locking Distal Lateral Femur Plate



Plate Thickness	6.0mm
Plate Shaft Width	15.40mm
Plate Head Width	30mm
Shaft Hole Spacing	20mm
Material	Titanium

Compatible with:
5.0mm Locking Screw
4.5 mm Cortical Screw
5.0mm Cancellous Screw
5.0mm Locking Cancellous Screw
5.0 mm Spacer

Ref. No	Hole	Direction	Plate Length
P00117	5	Left	156.00 mm
P00118	7	Left	196.00 mm
P00119	9	Left	236.00 mm
P00120	11	Left	276.00 mm
P00121	13	Left	316.00 mm
P00487	15	Left	356.00 mm
P00122	5	Right	156.00 mm
P00123	7	Right	196.00 mm
P00124	9	Right	236.00 mm
P00125	11	Right	276.00 mm
P00126	13	Right	316.00 mm
P00489	15	Right	356.00 mm

Locking Distal Medial Femur Plate



Material Titanium

Compatible with:
5.0mm Locking Screw
4.5 mm Cortical Screw
5.0mm Cancellous Screw
5.0mm Locking Cancellous Screw
5.0 mm Spacer

Ref. No	Hole	Direction
P00494	3	Left
P00495	5	Left
P00496	7	Left
P00497	9	Left
P00498	11	Left
P00499	3	Right
P00500	5	Right
P00501	7	Right
P00502	9	Right
P00503	11	Right

Locking Distal Femur (Lateral) Plate



Plate Thickness	4.0mm
Plate Shaft Width	13.60mm
Plate Head Width	22.50mm
Shaft Hole Spacing	
Material	Titanium

Compatible with:
5.0mm Locking Screw
4.5 mm Cortical Screw
5.0mm Cancellous Screw
5.0mm Locking Cancellous Screw
5.0 mm Spacer

Ref. No	Hole	Direction	Plate Length
P01190	4	Right	111.5 mm
P01191	4	Left	111.5 mm

4.5mm Locking Compression Plate Broad (LCP)



Plate Thickness	5.50mm
Plate Shaft Width	17.20mm
Plate Head Width	
Shaft Hole Spacing	18mm
Material	Titanium

Compatible with:
5.0mm Locking Screw
4.5 mm Cortical Screw
5.0mm Cancellous Screw
5.0mm Locking Cancellous Screw
5.0 mm Spacer

Ref. No	Hole	Plate Length
P00127	5	102 mm
P00128	6	120 mm
P00129	7	138 mm
P00130	8	156 mm
P00131	9	174 mm
P00132	10	192 mm
P00133	11	210 mm
P00134	12	228 mm
P00135	13	246 mm
P00136	14	264 mm
P00137	15	282 mm
P00138	16	300 mm

Knee

Patella Locking Plate




Plate Thickness	2.0mm
Plate Shaft Width	63mm
Plate Head Width	110mm
Shaft Hole Spacing	
Material	Titanium

Compatible with:
3.5mm Locking Screw
3.5mm Cortical Screw
4.0mm Cancellous Screw

Ref. No
P00801

HIP & PELVIC



Locking Proximal Femur Plate



Plate Thickness	6.0mm	Compatible with: 5.0 mm/6.5mm Locking Screw 4.5mm Cortical Screw 6.5mm Cancellous Screw 5.0 mm Spacer
Plate Shaft Width	17.50mm	
Plate Head Width	29mm	
Shaft Hole Spacing	18mm	
Material	Titanium	

Ref. No	Hole	Direction	Plate Length
P00159	4	Left	138.00 mm
P00160	6	Left	174.00 mm
P00161	8	Left	210.00 mm
P00162	10	Left	246.00 mm
P00163	12	Left	282.00 mm
P00504	14	Left	318.00mm
P00164	4	Right	138.00 mm
P00165	6	Right	174.00 mm
P00166	8	Right	210.00 mm
P00167	10	Right	246.00 mm
P00168	12	Right	282.00 mm
P00505	14	Right	318.00 mm

Locking 3.5 Pelvic Reconstruction Plate, Curved



Plate Thickness	3.0mm	Compatible with: 3.5mm Locking Screw 3.5mm Cortical Screw 4.0mm Cancellous Screw 4.0mm Locking Cancellous Screw 3.5 mm Spacer
Plate Shaft Width	10mm	
Plate Head Width		
Shaft Hole Spacing	Variable	
Material	Titanium	

Ref. No	Hole	Direction	Plate Length
P00591	4		60 mm
P00592	5		63.5 mm
P00258	6		72.4 mm
P00593	7		83.8 mm
P00259	8		95.1 mm
P00594	9		106 mm
P00260	10		116 mm
P00261	12		136 mm
P00262	14		154 mm
P00263	16		169.6 mm
P00595	18		183 mm
P00596	20		194 mm

Locking Reconstruction Plate




Plate Thickness	2.50mm	Compatible with: 3.5mm Locking Screw 3.5mm Cortical Screw 4.0mm Cancellous Screw 4.0mm Locking Cancellous Screw 3.5 mm Spacer
Plate Shaft Width	9.90mm	
Plate Head Width		
Shaft Hole Spacing	13mm	
Material	Titanium	

Ref. No	Hole	Plate Length
P00636	4	47 mm
P00637	5	60.6 mm
P00638	6	74.2 mm
P00639	7	87.8 mm
P00640	8	101.4 mm
P00641	9	115 mm
P00642	10	128.6 mm
P00643	12	155.8 mm
P00644	14	183 mm
P00645	16	210.2 mm
P00646	18	237.4 mm
P00647	20	264.6 mm

3.5mm Locking Reconstruction Plate, Curved



Plate Thickness	3.0mm	Compatible with: 3.5 mm Locking Screw 3.5 mm Cortical Screw 4.0mm Cancellous Screw 4.0mm Locking Cancellous Screw 3.5 mm Spacer
Plate Shaft Width	10mm	
Plate Head Width		
Shaft Hole Spacing	Variable	
Material	Titanium	

Ref. No	Hole	Plate Length
P00139	6	72.40 mm
P00140	7	83.80 mm
P00141	8	94.10 mm
P00142	9	106.00 mm
P00143	10	116.00 mm
P00144	12	136.00 mm
P00145	14	154.00 mm
P00146	16	169.60 mm
P00147	18	183.00 mm
P00148	20	194.00 mm

Locking Trochanteric Hook Plate



Plate Thickness	6.50mm
Plate Shaft Width	22mm
Plate Head Width	18.50mm
Shaft Hole Spacing	18mm
Material	Titanium

Compatible with:
Cable Tensioner
Cable Tensioner Screw
4.5 mm Cortical Screw
5.0 mm/6.5mm Locking Screw
6.5mm Cancellous Screw

Ref. No	Hole	Plate Length
P00651	2	45 mm
P00652	3	60 mm
P00653	4	130 mm
P00654	5	180 mm
P00655	6	213 mm
P00656	7	233 mm
P00657	8	265 mm
P00658	9	297 mm
P00659	10	329 mm

Locking Straight Trochanteric Plate



Material Titanium

Compatible with:
Cable Tensioner
Cable Tensioner Screw
4.5 mm Cortical Screw
5.0 mmLocking Screw

Ref. No	Hole
P01292	4
P01293	6
P01294	8

Button Implant



Material Titanium

Ref. No
S00800

Locking 5.0mm Connection Screw



Material Titanium

Ref. No
S00802

PEDIATRIC PLATES



3.5 mm Locking Curved Angled Head Plate with 3 Holes



Plate Thickness	3.0mm	Compatible with: 3.5mm Locking Screw 3.5mm Cortical Screw 4.0mm Cancellous Screw 4.0mm Locking Cancellous Screw 3.5 mm Spacer
Plate Shaft Width	12mm	
Plate Head Width	19.40mm	
Shaft Hole Spacing	13mm	
Material	Titanium	

Ref. No	Hole	Direction	Plate Length
P00610	3	100°	73 mm
P00611	3	110°	73 mm
P00612	3	120°	73 mm
P00613	3	140°	73 mm

Locking Curved Angled Head Plate with 3-Holes



Material	Titanium	Compatible with: 3.5mm Locking Screw 3.5mm Cortical Screw 4.0mm Cancellous Screw 4.0mm Locking Cancellous Screw 3.5 mm Spacer

Ref. No	Hole	Direction	Plate Length
P00614	3	100°	73 mm
P00615	3	110°	73 mm
P00616	3	120°	73 mm
P00617	3	140°	73 mm

4.5 mm Locking Curved Angled Head Plate with 3 Holes



Plate Thickness	3.50mm	Compatible with: 5.0mm Locking Screw 4.5 mm Cortical Screw 5.0mm Cancellous Screw 5.0mm Locking Cancellous Screw 5.0 mm Spacer
Plate Shaft Width	15.80mm	
Plate Head Width	23.30mm	
Shaft Hole Spacing		
Material	Titanium	

Ref. No	Hole	Direction	Plate Length
P00686	4	100°	106 mm
P00687	4	110°	106 mm
P00688	4	120°	106 mm
P00689	4	140°	106 mm

Eight Epiphyseal (Growth) Plate



Plate Thickness	1.50mm	Compatible with: 3.5mm Locking Screw 3.5mm Cortical Screw 4.0mm Cancellous Screw 4.5mm Cannulated Screw 3.5 mm Spacer
Plate Shaft Width	10mm	
Plate Head Width	10mm	
Shaft Hole Spacing	2mm	
Material	Titanium	

Ref. No	Hole	Plate Length
P00713	12 mm	28 mm
P00714	16 mm	32 mm
P00715	18 mm	34 mm

SCREWS & ACCESSORIES

Cortex Screw, Self-Tapping Ø 3.5 mm



Material Titanium



Locking Screw, Self-Tapping Ø 3.5 mm



Material Titanium




Ref. No	Diameter	Length
S00955	Ø 3.5 mm	10 mm
S00140	Ø 3.5 mm	12 mm
S00141	Ø 3.5 mm	14 mm
S00142	Ø 3.5 mm	16 mm
S00143	Ø 3.5 mm	18 mm
S00144	Ø 3.5 mm	20 mm
S00145	Ø 3.5 mm	22 mm
S00146	Ø 3.5 mm	24 mm
S00147	Ø 3.5 mm	26 mm
S00148	Ø 3.5 mm	28 mm
S00149	Ø 3.5 mm	30 mm
S00150	Ø 3.5 mm	32 mm
S00151	Ø 3.5 mm	34 mm
S00152	Ø 3.5 mm	36 mm
S00153	Ø 3.5 mm	38 mm
S00154	Ø 3.5 mm	40 mm
S00155	Ø 3.5 mm	42 mm
S00156	Ø 3.5 mm	44 mm
S00157	Ø 3.5 mm	46 mm
S00158	Ø 3.5 mm	48 mm
S00159	Ø 3.5 mm	50 mm
S00559	Ø 3.5 mm	55 mm
S00560	Ø 3.5 mm	60 mm
S00561	Ø 3.5 mm	65 mm
S00562	Ø 3.5 mm	70 mm
S00563	Ø 3.5 mm	75 mm
S00564	Ø 3.5 mm	80 mm
S00565	Ø 3.5 mm	85 mm
S00566	Ø 3.5 mm	90 mm
S00567	Ø 3.5 mm	95 mm
S00568	Ø 3.5 mm	100 mm
S00569	Ø 3.5 mm	105 mm
S00570	Ø 3.5 mm	110 mm
S00956	Ø 3.5 mm	115 mm
S00571	Ø 3.5 mm	120 mm
S00572	Ø 3.5 mm	125 mm
S00573	Ø 3.5 mm	130 mm
S00574	Ø 3.5 mm	135 mm
S00575	Ø 3.5 mm	140 mm

Ref. No	Diameter	Length
S00196	Ø 3.5 mm	10 mm
S00197	Ø 3.5 mm	12 mm
S00198	Ø 3.5 mm	14 mm
S00199	Ø 3.5 mm	16 mm
S00200	Ø 3.5 mm	18 mm
S00201	Ø 3.5 mm	20 mm
S00202	Ø 3.5 mm	22 mm
S00203	Ø 3.5 mm	24 mm
S00204	Ø 3.5 mm	26 mm
S00205	Ø 3.5 mm	28 mm
S00206	Ø 3.5 mm	30 mm
S00207	Ø 3.5 mm	32 mm
S00208	Ø 3.5 mm	34 mm
S00209	Ø 3.5 mm	35 mm
S00210	Ø 3.5 mm	36 mm
S00211	Ø 3.5 mm	38 mm
S00212	Ø 3.5 mm	40 mm
S00213	Ø 3.5 mm	42 mm
S00214	Ø 3.5 mm	44 mm
S00215	Ø 3.5 mm	46 mm
S00216	Ø 3.5 mm	48 mm
S00217	Ø 3.5 mm	50 mm
S00218	Ø 3.5 mm	55 mm
S00219	Ø 3.5 mm	60 mm
S00220	Ø 3.5 mm	65 mm
S00577	Ø 3.5 mm	70 mm
S00578	Ø 3.5 mm	75 mm
S00579	Ø 3.5 mm	80 mm
S00580	Ø 3.5 mm	85 mm
S00581	Ø 3.5 mm	90 mm
S00582	Ø 3.5 mm	95 mm
S00583	Ø 3.5 mm	100 mm
S00584	Ø 3.5 mm	105 mm
S00585	Ø 3.5 mm	110 mm
S00586	Ø 3.5 mm	115 mm
S00587	Ø 3.5 mm	120 mm
S00588	Ø 3.5 mm	125 mm
S00589	Ø 3.5 mm	130 mm
S00590	Ø 3.5 mm	135 mm
S00591	Ø 3.5 mm	140 mm

Cortex Screw Ø 4.5 mm


Material Titanium



Locking Screw, Self-Tapping Ø 5.0 mm


Material Titanium


Ref. No	Diameter	Length
S00592	Ø 4.5 mm	14 mm
S00290	Ø 4.5 mm	16 mm
S00291	Ø 4.5 mm	18 mm
S00292	Ø 4.5 mm	20 mm
S00293	Ø 4.5 mm	22 mm
S00294	Ø 4.5 mm	24 mm
S00295	Ø 4.5 mm	26 mm
S00296	Ø 4.5 mm	28 mm
S00297	Ø 4.5 mm	30 mm
S00298	Ø 4.5 mm	32 mm
S00299	Ø 4.5 mm	34 mm
S00300	Ø 4.5 mm	36 mm
S00301	Ø 4.5 mm	38 mm
S00302	Ø 4.5 mm	40 mm
S00303	Ø 4.5 mm	42 mm
S00304	Ø 4.5 mm	44 mm
S00305	Ø 4.5 mm	46 mm
S00306	Ø 4.5 mm	48 mm
S00307	Ø 4.5 mm	50 mm
S00308	Ø 4.5 mm	52 mm
S00309	Ø 4.5 mm	54 mm
S00310	Ø 4.5 mm	56 mm
S00311	Ø 4.5 mm	58 mm
S00312	Ø 4.5 mm	60 mm
S00593	Ø 4.5 mm	65 mm
S00594	Ø 4.5 mm	70 mm
S00595	Ø 4.5 mm	75 mm
S00596	Ø 4.5 mm	80 mm
S00597	Ø 4.5 mm	85 mm
S00598	Ø 4.5 mm	90 mm
S00599	Ø 4.5 mm	95 mm
S00600	Ø 4.5 mm	100 mm
S00601	Ø 4.5 mm	105 mm
S00602	Ø 4.5 mm	110 mm
S00603	Ø 4.5 mm	115 mm
S00604	Ø 4.5 mm	120 mm
S00605	Ø 4.5 mm	125 mm
S00606	Ø 4.5 mm	130 mm
S00607	Ø 4.5 mm	135 mm
S00608	Ø 4.5 mm	140 mm

Ref. No	Diameter	Length
S00313	Ø 5.0 mm	14 mm
S00314	Ø 5.0 mm	16 mm
S00315	Ø 5.0 mm	18 mm
S00316	Ø 5.0 mm	20 mm
S00317	Ø 5.0 mm	22 mm
S00318	Ø 5.0 mm	24 mm
S00319	Ø 5.0 mm	26 mm
S00320	Ø 5.0 mm	28 mm
S00321	Ø 5.0 mm	30 mm
S00322	Ø 5.0 mm	32 mm
S00323	Ø 5.0 mm	34 mm
S00324	Ø 5.0 mm	36 mm
S00325	Ø 5.0 mm	38 mm
S00326	Ø 5.0 mm	40 mm
S00327	Ø 5.0 mm	42 mm
S00328	Ø 5.0 mm	44 mm
S00329	Ø 5.0 mm	45 mm
S00330	Ø 5.0 mm	46 mm
S00331	Ø 5.0 mm	48 mm
S00332	Ø 5.0 mm	50 mm
S00333	Ø 5.0 mm	55 mm
S00334	Ø 5.0 mm	60 mm
S00335	Ø 5.0 mm	65 mm
S00336	Ø 5.0 mm	70 mm
S00337	Ø 5.0 mm	75 mm
S00338	Ø 5.0 mm	80 mm
S00339	Ø 5.0 mm	85 mm
S00340	Ø 5.0 mm	90 mm
S00341	Ø 5.0 mm	95 mm
S00342	Ø 5.0 mm	100 mm
S00609	Ø 5.0 mm	105 mm
S00610	Ø 5.0 mm	110mm
S00611	Ø 5.0 mm	115mm
S00612	Ø 5.0 mm	120mm


Locking Cancellous Screw, Self-Tapping Ø 6.5 mm



Material Titanium



Cannulated Locking Screw, Self-Tapping Ø 6.5 mm



Material Titanium



Ref. No	Diameter	Length
S00748	Ø 6.5 mm	30 mm
S00749	Ø 6.5 mm	35 mm
S00750	Ø 6.5 mm	40 mm
S00751	Ø 6.5 mm	45 mm
S00752	Ø 6.5 mm	50 mm
S00753	Ø 6.5 mm	55 mm
S00754	Ø 6.5 mm	60 mm
S00755	Ø 6.5 mm	65 mm
S00756	Ø 6.5 mm	70 mm
S00757	Ø 6.5 mm	75 mm
S00758	Ø 6.5 mm	80 mm
S00759	Ø 6.5 mm	85 mm
S00760	Ø 6.5 mm	90 mm
S00761	Ø 6.5 mm	95 mm
S00762	Ø 6.5 mm	100 mm
S00763	Ø 6.5 mm	105 mm
S00764	Ø 6.5 mm	110 mm
S00765	Ø 6.5 mm	115 mm
S00766	Ø 6.5 mm	120 mm

Ref. No	Diameter	Length
S00522	Ø 6.5 mm	50 mm
S00523	Ø 6.5 mm	55 mm
S00524	Ø 6.5 mm	60 mm
S00525	Ø 6.5 mm	65 mm
S00526	Ø 6.5 mm	70 mm
S00527	Ø 6.5 mm	75 mm
S00528	Ø 6.5 mm	80 mm
S00529	Ø 6.5 mm	85 mm
S00530	Ø 6.5 mm	90 mm
S00531	Ø 6.5 mm	95 mm
S00532	Ø 6.5 mm	100 mm
S00533	Ø 6.5 mm	105 mm
S00534	Ø 6.5 mm	110 mm
S00535	Ø 6.5 mm	115 mm
S00536	Ø 6.5 mm	120 mm

Cancellous Bone Screw, Partially Threaded Ø 4.0 mm



Material Titanium



Ref. No	Diameter	Length
S00160	Ø 4.0 mm	20 mm
S00161	Ø 4.0 mm	22 mm
S00162	Ø 4.0 mm	24 mm
S00163	Ø 4.0 mm	26 mm
S00164	Ø 4.0 mm	28 mm
S00165	Ø 4.0 mm	30 mm
S00166	Ø 4.0 mm	32 mm
S00167	Ø 4.0 mm	34 mm
S00625	Ø 4.0 mm	35 mm
S00168	Ø 4.0 mm	36 mm
S00169	Ø 4.0 mm	38 mm
S00170	Ø 4.0 mm	40 mm
S00171	Ø 4.0 mm	42 mm
S00172	Ø 4.0 mm	44 mm
S00626	Ø 4.0 mm	45 mm
S00173	Ø 4.0 mm	46 mm
S00174	Ø 4.0 mm	48 mm
S00175	Ø 4.0 mm	50 mm
S00627	Ø 4.0 mm	55 mm
S00628	Ø 4.0 mm	60 mm
S00629	Ø 4.0 mm	65 mm
S00630	Ø 4.0 mm	70 mm

Cancellous Bone Screw, Fully Threaded Ø 4.0 mm



Material Titanium



Ref. No	Diameter	Length
S00176	Ø 4.0 mm	12 mm
S00177	Ø 4.0 mm	14 mm
S00178	Ø 4.0 mm	16 mm
S00179	Ø 4.0 mm	18 mm
S00180	Ø 4.0 mm	20 mm
S00181	Ø 4.0 mm	22 mm
S00182	Ø 4.0 mm	24 mm
S00183	Ø 4.0 mm	26 mm
S00184	Ø 4.0 mm	28 mm
S00185	Ø 4.0 mm	30 mm
S00186	Ø 4.0 mm	32 mm
S00187	Ø 4.0 mm	34 mm
S00631	Ø 4.0 mm	35 mm
S00188	Ø 4.0 mm	36 mm
S00189	Ø 4.0 mm	38 mm
S00190	Ø 4.0 mm	40 mm
S00191	Ø 4.0 mm	42 mm
S00192	Ø 4.0 mm	44 mm
S00632	Ø 4.0 mm	45 mm
S00193	Ø 4.0 mm	46 mm
S00194	Ø 4.0 mm	48 mm
S00195	Ø 4.0 mm	50 mm
S00633	Ø 4.0 mm	55 mm
S00634	Ø 4.0 mm	60 mm
S00635	Ø 4.0 mm	65 mm
S00636	Ø 4.0 mm	70 mm

Cancellous Screw Partially Threaded Ø 5.0 mm



Material Titanium



Ref. No	Diameter	Length
S00459	Ø 5.0 mm	14 mm
S00460	Ø 5.0 mm	16 mm
S00461	Ø 5.0 mm	18 mm
S00462	Ø 5.0 mm	20 mm
S00463	Ø 5.0 mm	22 mm
S00464	Ø 5.0 mm	24 mm
S00465	Ø 5.0 mm	26 mm
S00466	Ø 5.0 mm	28 mm
S00467	Ø 5.0 mm	30 mm
S00468	Ø 5.0 mm	32 mm
S00469	Ø 5.0 mm	34 mm
S00470	Ø 5.0 mm	36 mm
S00471	Ø 5.0 mm	38 mm
S00472	Ø 5.0 mm	40 mm
S00473	Ø 5.0 mm	42 mm
S00474	Ø 5.0 mm	44 mm
S00475	Ø 5.0 mm	45 mm

Ref. No	Diameter	Length
S00476	Ø 5.0 mm	46 mm
S00477	Ø 5.0 mm	48 mm
S00478	Ø 5.0 mm	50 mm
S00479	Ø 5.0 mm	55 mm
S00480	Ø 5.0 mm	60 mm
S00481	Ø 5.0 mm	65 mm
S00482	Ø 5.0 mm	70 mm
S00483	Ø 5.0 mm	75 mm
S00484	Ø 5.0 mm	80 mm
S00485	Ø 5.0 mm	85 mm
S00486	Ø 5.0 mm	90 mm
S00487	Ø 5.0 mm	95 mm
S00488	Ø 5.0 mm	100 mm
S00637	Ø 5.0 mm	105 mm
S00638	Ø 5.0 mm	110 mm
S00639	Ø 5.0 mm	115 mm
S00640	Ø 5.0 mm	120 mm

Cancellous Screw Fully Threaded 5.0mm



Material Titanium



Ref. No	Diameter	Length
S00839	Ø 5.0 mm	14mm
S00840	Ø 5.0 mm	16mm
S00841	Ø 5.0 mm	18mm
S00842	Ø 5.0 mm	20mm
S00843	Ø 5.0 mm	22mm
S00844	Ø 5.0 mm	24mm
S00845	Ø 5.0 mm	26mm
S00846	Ø 5.0 mm	28mm
S00847	Ø 5.0 mm	30mm
S00848	Ø 5.0 mm	32mm
S00849	Ø 5.0 mm	34mm
S00850	Ø 5.0 mm	36mm
S00851	Ø 5.0 mm	38mm
S00852	Ø 5.0 mm	40mm
S00853	Ø 5.0 mm	42mm
S00854	Ø 5.0 mm	44mm
S00855	Ø 5.0 mm	45mm
S00856	Ø 5.0 mm	46mm
S00857	Ø 5.0 mm	48mm
S00858	Ø 5.0 mm	50mm
S00859	Ø 5.0 mm	55mm
S00860	Ø 5.0 mm	60mm
S00861	Ø 5.0 mm	65mm
S00862	Ø 5.0 mm	70mm
S00863	Ø 5.0 mm	75mm
S00864	Ø 5.0 mm	80mm
S00865	Ø 5.0 mm	85mm
S00866	Ø 5.0 mm	90mm
S00867	Ø 5.0 mm	95mm
S00868	Ø 5.0 mm	100mm
S00869	Ø 5.0 mm	105mm
S00870	Ø 5.0 mm	110mm
S00871	Ø 5.0 mm	115mm
S00872	Ø 5.0 mm	120mm

Cancellous Bone Screw, Half Threaded Ø 6.5 mm



Material Titanium



Ref. No	Diameter	Length
S00645	Ø 6.5 mm	30 mm
S00646	Ø 6.5 mm	35 mm
S00647	Ø 6.5 mm	40 mm
S00648	Ø 6.5 mm	45 mm
S00649	Ø 6.5 mm	50 mm
S00650	Ø 6.5 mm	55 mm
S00651	Ø 6.5 mm	60 mm
S00652	Ø 6.5 mm	65 mm
S00653	Ø 6.5 mm	70 mm
S00654	Ø 6.5 mm	75 mm
S00655	Ø 6.5 mm	80 mm
S00656	Ø 6.5 mm	85 mm
S00657	Ø 6.5 mm	90 mm
S00658	Ø 6.5 mm	95 mm
S00659	Ø 6.5 mm	100 mm
S00660	Ø 6.5 mm	105 mm
S00661	Ø 6.5 mm	110 mm
S00662	Ø 6.5 mm	115 mm
S00663	Ø 6.5 mm	120 mm

Cancellous Bone Screw, Fully Threaded Ø 6.5 mm



Material Titanium



Ref. No	Diameter	Length
S00641	Ø 6.5 mm	30 mm
S00642	Ø 6.5 mm	35 mm
S00643	Ø 6.5 mm	40 mm
S00644	Ø 6.5 mm	45 mm
S00489	Ø 6.5 mm	50 mm
S00490	Ø 6.5 mm	55 mm
S00491	Ø 6.5 mm	60 mm
S00492	Ø 6.5 mm	65 mm
S00493	Ø 6.5 mm	70 mm
S00494	Ø 6.5 mm	75 mm
S00495	Ø 6.5 mm	80 mm
S00496	Ø 6.5 mm	85 mm
S00497	Ø 6.5 mm	90 mm
S00498	Ø 6.5 mm	95 mm
S00499	Ø 6.5 mm	100 mm
S00500	Ø 6.5 mm	105 mm
S00501	Ø 6.5 mm	110 mm
S00502	Ø 6.5 mm	115 mm
S00503	Ø 6.5 mm	120 mm


4.0mm Locking, Cancellous Screw, Fully Threaded



Material Titanium



Locking Cancellous Screw, Self-Tapping Ø 5.0 mm




Material Titanium



Ref. No	Diameter	Length
S00684	Ø 4.0 mm	10 mm
S00685	Ø 4.0 mm	12 mm
S00686	Ø 4.0 mm	14 mm
S00687	Ø 4.0 mm	16 mm
S00688	Ø 4.0 mm	18 mm
S00689	Ø 4.0 mm	20 mm
S00690	Ø 4.0 mm	22 mm
S00691	Ø 4.0 mm	24 mm
S00692	Ø 4.0 mm	26 mm
S00693	Ø 4.0 mm	28 mm
S00694	Ø 4.0 mm	30 mm
S00695	Ø 4.0 mm	35 mm
S00696	Ø 4.0 mm	40 mm
S00697	Ø 4.0 mm	45 mm
S00698	Ø 4.0 mm	50 mm
S00699	Ø 4.0 mm	55 mm
S00700	Ø 4.0 mm	60 mm
S00701	Ø 4.0 mm	65 mm
S00702	Ø 4.0 mm	70 mm
S00703	Ø 4.0 mm	75 mm
S00704	Ø 4.0 mm	80 mm
S00705	Ø 4.0 mm	85 mm
S00706	Ø 4.0 mm	90 mm
S00707	Ø 4.0 mm	95 mm
S00708	Ø 4.0 mm	100 mm
S00709	Ø 4.0 mm	105 mm
S00710	Ø 4.0 mm	110 mm
S00711	Ø 4.0 mm	115 mm
S00712	Ø 4.0 mm	120 mm
S00713	Ø 4.0 mm	125 mm
S00714	Ø 4.0 mm	130 mm
S00715	Ø 4.0 mm	135 mm
S00716	Ø 4.0 mm	140 mm

Ref. No	Diameter	Length
S00429	Ø 5.0 mm	14 mm
S00430	Ø 5.0 mm	16 mm
S00431	Ø 5.0 mm	18 mm
S00432	Ø 5.0 mm	20 mm
S00433	Ø 5.0 mm	22 mm
S00434	Ø 5.0 mm	24 mm
S00435	Ø 5.0 mm	26 mm
S00436	Ø 5.0 mm	28 mm
S00437	Ø 5.0 mm	30 mm
S00438	Ø 5.0 mm	32 mm
S00439	Ø 5.0 mm	34 mm
S00440	Ø 5.0 mm	36 mm
S00441	Ø 5.0 mm	38 mm
S00442	Ø 5.0 mm	40 mm
S00443	Ø 5.0 mm	42 mm
S00444	Ø 5.0 mm	44 mm
S00445	Ø 5.0 mm	45 mm
S00446	Ø 5.0 mm	46 mm
S00447	Ø 5.0 mm	48 mm
S00448	Ø 5.0 mm	50 mm
S00449	Ø 5.0 mm	55 mm
S00450	Ø 5.0 mm	60 mm
S00451	Ø 5.0 mm	65 mm
S00452	Ø 5.0 mm	70 mm
S00453	Ø 5.0 mm	75 mm
S00454	Ø 5.0 mm	80 mm
S00455	Ø 5.0 mm	85 mm
S00456	Ø 5.0 mm	90 mm
S00457	Ø 5.0 mm	95 mm
S00458	Ø 5.0 mm	100 mm
S00740	Ø 5.0 mm	105 mm
S00741	Ø 5.0 mm	110 mm
S00742	Ø 5.0 mm	115 mm
S00743	Ø 5.0 mm	120 mm
S00744	Ø 5.0 mm	125 mm
S00745	Ø 5.0 mm	130 mm
S00746	Ø 5.0 mm	135 mm
S00747	Ø 5.0 mm	140 mm

5.0mm Locking Spacer Screw



Material Titanium

Ref. No	Diameter
S00799	5.0mm

3.5mm Locking Spacer Screw



Material Titanium



Ref. No	Diameter
S00798	3.5mm

Washer



Material Titanium

Ref. No	
PL00001	Small
PL00002	Large

A black and white photograph of a tall, modern glass skyscraper with a distinctive stepped-top design. The building is surrounded by dark, dense foliage at the bottom. A red diagonal line runs across the lower portion of the image, separating the text area from the bottom section.

We blend manufacturing flexibility with a deep understanding of the global medical market. We know that in surgery, reliability is everything. That is why we are dedicated to delivering high-performance implants for Orthopedic Trauma, Spine, and Craniomaxillofacial (CMF) surgery that you can trust case after case.

Our goal is simple; to support healthcare providers and distributors with products that offer clinical excellence, reliable supply, and genuine economic value.

We have built our product line to meet the realities of modern surgical practice. Our core areas of expertise include:

Orthopedic Trauma - A complete range of locking plates, cannulated screws, and intramedullary nailing systems.

Spine - Comprehensive solutions including stabilization systems, PEEK cages, cervical plates, and disc prostheses.

Craniomaxillofacial (CMF) - Specialized plating and screw systems designed for precision in cranial and facial reconstruction.

Quality and consistency are the heart of our operation. Oltho operates under a Quality Management System certified to ISO 13485:2016. By using the latest production technology and strict control processes, we ensure every implant meets international regulatory standards and CE compliance. We ensure regulatory continuity and supply stability by maintaining our CE certification in accordance with the EU MDR extension.





oltho.com

info@oltho.com

Tel . +90 312 441 6265

Fax. +90 312 441 6267

Kızılırmak Mahallesi
Muhsin Yazıcıoğlu Caddesi
No: 39/A-57 Çankaya
Ankara - **TÜRKİYE**

oltho

oltho
ITALIANO

TRAUMA

For our complete range of systems, please tab or scan the QR code:

