



UTC

Drug Tester

Testing 7 drug types

Easy data entry and menu selection using keypad built in the device

Capability of saving driver's name-surname, driver's ID number, Police's ID number, police's registration number, police's name-surname in its internal memory

Capability to display the status of the memory and search for records

Capability to transfer recorded data from device to computer via USB with UTC-DATA software

Printing reports from portable printer via Bluetooth communication

Record/report of refusal feature for those who refuse to take a measurement

Recorded data in the device cannot be deleted by the user

Ambient temperature control

Restriction of the personal authorization

Possibility of easy to use due to the wide graphic display

Rechargeable Li-ion battery

Built-in battery for clock and date

220VAC and 12VDC charging kits

CE certification, compliance with EMC and low voltage directives


Keypad illumination and sound on/off

Display contrast adjustment



UTC-DATA Software for Data Transfer;

With UTC-DATA software, you can transfer data from device to computer, archive, search records, transfer records to different media (memory cards, floppy disks, etc.) and printing



Carrying case



UTK

Drug Tester

Easy-to-use saliva collection pad integrated test kit

Control lines indicating that the sufficient amount of sample is collected on the test kit

Marking on the package indicating lot, date of manufacture / expiry date and type of drugs

Built-in saliva collection pad

Optional RF ID labeling

Compliance with EMC and low voltage directive

DRUG KINDS and THRESHOLD VALUES

DRUGS	THRESHOLD (ng/ml)
Amphetamine (AMP) Methamphetamine (MET) Ecstasy (MDMA)	50
Benzodiazepines (BZD)	20
Cocaine (COC)	30
Opiates (OPI)	40
Marijuana (THC)	25



PRN-08 Printer

Printing method	Direct Thermal
Battery	Rechargeable Li-ion or Li-polymer
Interface	Bluetooth
Paper width	57,5 mm
Diameter of paper roll	40 mm
Indicators	LED (error, no paper, open cover, charge, normal/low battery)
Automatic shutdown	Yes
Self-Test	Yes

